



**November 8, 2024**

**Nevada Medicaid Web Announcement 3477**

## **Insulin Pump Recall: Medtronic Notifies Users of MiniMed 600 and 700 Series Pumps of Risk of Shorter than Expected Battery Life**

Medtronic is notifying users of MiniMed 600 and 700 series insulin pumps of an increased risk for reduced battery life and less time until shutdown after a battery alert occurs. Pumps that have been dropped, bumped, or experienced another physical impact may have damaged electrical components causing this issue—even a single drop can impact battery life.

The U.S. Food and Drug Administration (FDA) has identified this as a Class I recall, the most serious type of recall. This recall does not involve removing all products. Users should contact Medtronic if they experience significant reduced battery life to determine if a new pump is needed.

See the [full recall notice](#) on the FDA's website for more information and the next actions to take.

Customers in the U.S. can call Medtronic at 1-800-378-2292 for 24-hour technical support.